

PPM1J Polyclonal Antibody

Cat No: HR1AP11170

For research use only

Overview

Product Name	PPM1J Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	PPM1J PPP2CZ
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	protein phosphatase, Mg ²⁺ /Mn ²⁺ dependent 1J(PPM1J) Homo sapiens This gene encodes the serine/threonine protein phosphatase. The mouse homolog of this gene apparently belongs to the protein phosphatase 2C family of genes. The exact function of this gene is not yet known. [provided by RefSeq, Jul 2008],