

# BD1L1 Polyclonal Antibody

Cat No: HR1AP11165

For research use only

## Overview

Product Name	BD1L1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	335kDa
GeneID?Human?	BOD1L1 BOD1L FAM44A KIAA1327
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Component of the fork protection machinery required to protect stalled/damaged replication forks from uncontrolled DNA2-dependent resection. Acts by stabilizing RAD51 at stalled replication forks and protecting RAD51 nucleofilaments from the antirecombinogenic activities of FBXO18/FBH1 and BLM.