

# FOH1B Polyclonal Antibody

Cat No: HR1AP11164

For research use only

## Overview

Product Name	FOH1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	FOLH1B PSMAL GIG26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	FOLH1B (Folate Hydrolase 1B) is a Protein Coding gene. Among its related pathways are Metabolism and Amino acid synthesis and interconversion (transamination). Has both folate hydrolase and N-acetylated-alpha-linked-acidic dipeptidase (NAALADase) activity. Exhibits a dipeptidyl-peptidase IV type activity.