

# POTEE Polyclonal Antibody

Cat No: HR1AP11163

For research use only

## Overview

Product Name	POTEE Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	118kDa
GeneID?Human?	POTEE A26C1A POTE2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the actin family.,similarity:Contains 5 ANK repeats.,similarity:In the C-terminal section; belongs to the actin family.,similarity:In the N-terminal section; belongs to the POTE family.,