

## PNKD Polyclonal Antibody

Cat No: HR1AP11161

For research use only

### Overview

Product Name	PNKD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	PNKD KIAA1184 MR1 TAHCCP2 FKSG19 UNQ2491/PRO5778
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	paroxysmal nonkinesigenic dyskinesia(PNKD) Homo sapiens This gene is thought to play a role in the regulation of myofibrillogenesis. Mutations in this gene have been associated with the movement disorder paroxysmal non-kinesigenic dyskinesia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010],