

## SYF2 Polyclonal Antibody

Cat No: HR1AP11158

For research use only

### Overview

Product Name	SYF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	SYF2 CBPIN GCIPIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SYF2 pre-mRNA splicing factor(SYF2) Homo sapiens This gene encodes a nuclear protein that interacts with cyclin D-type binding-protein 1, which is thought to be a cell cycle regulator at the G1/S transition. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],