

PRP8 Polyclonal Antibody

Cat No: HR1AP11155

For research use only

Overview

Product Name	PRP8 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	256kDa
GeneID?Human?	PRPF8 PRPC8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>pre-mRNA processing factor 8(PRPF8) Homo sapiens Pre-mRNA splicing occurs in 2 sequential transesterification steps. The protein encoded by this gene is a component of both U2- and U12-dependent spliceosomes, and found to be essential for the catalytic step II in pre-mRNA splicing process. It contains several WD repeats, which function in protein-protein interactions. This protein has a sequence similarity to yeast Prp8 protein. This gene is a candidate gene for autosomal dominant retinitis pigmentosa. [provided by RefSeq, Jul 2008].</p>