

PR40A Polyclonal Antibody

Cat No: HR1AP11154

For research use only

Overview

Product Name	PR40A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	PRPF40A FBP11 FLAF1 FBNP3 HIP10 HYPA HSPC225
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>domain:The WW domains are essential for localization to nuclear speckles.,function:Binds to WASL/N-WASP and suppresses its translocation from the nucleus to the cytoplasm, thereby inhibiting its cytoplasmic function (By similarity). May be involved in pre-mRNA splicing.,sequence caution:Contaminating sequence. Potential poly-A sequence starting in position 409.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Belongs to the PRPF40 family.,similarity:Contains 2 WW domains.,similarity:Contains 5 FF domains.,subunit:Interacts with the N-terminus of HTT. Interacts with the phosphorylated carboxyterminal domain of POLR2A. Interacts with SF1, SRPK1 CARD8, ATBF1 and MECP2 (By similarity). Interacts through the WW domains with formin proline-rich regions and with WASL/N-WASP.,tissue specificity:Widely expressed.,</p>