

# PCF11 Polyclonal Antibody

Cat No: HR1AP11152

For research use only

## Overview

Product Name	PCF11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	171kDa
GeneID?Human?	PCF11 KIAA0824
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PCF11 cleavage and polyadenylation factor subunit(PCF11) Homo sapiens The protein encoded by this gene binds to CLP1 to form pre-mRNA cleavage factor IIm. The encoded protein is necessary for efficient Pol II transcription termination and may be involved in degradation of the 3' product of polyA site cleavage. [provided by RefSeq, Oct 2016],