

## HCN3 Polyclonal Antibody

Cat No: HR1AP11147

For research use only

### Overview

Product Name	HCN3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	85kDa
GeneID?Human?	HCN3 KIAA1535
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hyperpolarization activated cyclic nucleotide gated potassium channel 3(HCN3) Homo sapiens This gene encodes a multi-pass membrane protein that functions as a voltage gated cation channel. The encoded protein is a member of a family of closely related cyclic adenosine monophosphate-binding channel proteins. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012],