

## HCN1 Polyclonal Antibody

Cat No: HR1AP11146

For research use only

### Overview

Product Name	HCN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	97kDa
GeneID?Human?	HCN1 BCNG1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>hyperpolarization activated cyclic nucleotide gated potassium channel 1(HCN1) Homo sapiens The membrane protein encoded by this gene is a hyperpolarization-activated cation channel that contributes to the native pacemaker currents in heart and neurons. The encoded protein can homodimerize or heterodimerize with other pore-forming subunits to form a potassium channel. This channel may act as a receptor for sour tastes. [provided by RefSeq, Oct 2011].</p>