

PAIP1 Polyclonal Antibody

Cat No: HR1AP11129

For research use only

Overview

Product Name	PAIP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	PAIP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>poly(A) binding protein interacting protein 1(PAIP1) Homo sapiens The protein encoded by this gene interacts with poly(A)-binding protein and with the cap-binding complex eIF4A. It is involved in translational initiation and protein biosynthesis. Overexpression of this gene in COS7 cells stimulates translation. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. [provided by RefSeq, Jul 2008],</p>