

PLXA2 Polyclonal Antibody

Cat No: HR1AP11124

For research use only

Overview

Product Name	PLXA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	208kDa
GeneID?Human?	PLXNA2 KIAA0463 OCT PLXN2 UNQ209/PRO235
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>plexin A2(PLXNA2) Homo sapiens This gene encodes a member of the plexin-A family of semaphorin co-receptors. Semaphorins are a large family of secreted or membrane-bound proteins that mediate repulsive effects on axon pathfinding during nervous system development. A subset of semaphorins are recognized by plexin-A/neuropilin transmembrane receptor complexes, triggering a cellular signal transduction cascade that leads to axon repulsion. This plexin-A family member is thought to transduce signals from semaphorin-3A and -3C. [provided by RefSeq, Jul 2008],</p>