

PAIRB Polyclonal Antibody

Cat No: HR1AP1121

For research use only

Overview

Product Name	PAIRB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	SERBP1 PAIRBP1 CGI-55
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in the regulation of mRNA stability. Binds to the 3'-most 134 nt of the SERPINE1/PAI1 mRNA, a region which confers cyclic nucleotide regulation of message decay.,subcellular location:Also found in perinuclear regions.,subunit:Interacts with CHD3 and TDRD3.,tissue specificity:Expressed at high level in the heart, skeletal muscle and kidney, and at low levels in placenta, liver and brain.,