

# PKHG1 Polyclonal Antibody

Cat No: HR1AP11122

For research use only

## Overview

Product Name	PKHG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	152kDa
GeneID?Human?	PLEKHG1 KIAA1209
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,