

## PIGW Polyclonal Antibody

Cat No: HR1AP11107

For research use only

### Overview

Product Name	PIGW Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	PIGW
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>phosphatidylinositol glycan anchor biosynthesis class W(PIGW) Homo sapiens Glycosylphosphatidylinositol (GPI) is a complex glycolipid that anchors many proteins to the cell surface. PIGW acts in the third step of GPI biosynthesis and acylates the inositol ring of phosphatidylinositol (Murakami et al., 2003 [PubMed 14517336]).[supplied by OMIM, Mar 2008],</p>