

PI5L1 Polyclonal Antibody

Cat No: HR1AP11105

For research use only

Overview

Product Name	PI5L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	PIP5KL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	phosphatidylinositol-4-phosphate 5-kinase like 1(PIP5KL1) Homo sapiens PIP5KL1 is a phosphoinositide kinase-like protein that lacks intrinsic lipid kinase activity but associates with type I PIPKs (see PIP5K1A; MIM 603275) and may play a role in localization of PIPK activity (Chang et al., 2004 [PubMed 14701839]).[supplied by OMIM, Jun 2009],