

TPTE2 Polyclonal Antibody

Cat No: HR1AP11104

For research use only

Overview

Product Name	TPTE2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	TPTE2 TPIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transmembrane phosphoinositide 3-phosphatase and tensin homolog 2(TPTE2) Homo sapiens TPIP is a member of a large class of membrane-associated phosphatases with substrate specificity for the 3-position phosphate of inositol phospholipids.[supplied by OMIM, Jul 2002],