

PREX1 Polyclonal Antibody

Cat No: HR1AP11098

For research use only

Overview

Product Name	PREX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	182kDa
GeneID?Human?	PREX1 KIAA1415
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>phosphatidylinositol-3,4,5-trisphosphate dependent Rac exchange factor 1(PREX1) Homo sapiens The protein encoded by this gene acts as a guanine nucleotide exchange factor for the RHO family of small GTP-binding proteins (RACs). It has been shown to bind to and activate RAC1 by exchanging bound GDP for free GTP. The encoded protein, which is found mainly in the cytoplasm, is activated by phosphatidylinositol-3,4,5-trisphosphate and the beta-gamma subunits of heterotrimeric G proteins. [provided by RefSeq, Jul 2008],</p>