

PM34 Polyclonal Antibody

Cat No: HR1AP11092

For research use only

Overview

Product Name	PM34 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	SLC25A17 PMP34
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 25 member 17(SLC25A17) Homo sapiens This gene encodes a peroxisomal membrane protein that belongs to the family of mitochondrial solute carriers. It is expressed in the liver, and is likely involved in transport. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2013],