

PEX13 Polyclonal Antibody

Cat No: HR1AP11090

For research use only

Overview

Product Name	PEX13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	PEX13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>peroxisomal biogenesis factor 13(PEX13) Homo sapiens This gene encodes a peroxisomal membrane protein that binds the type 1 peroxisomal targeting signal receptor via a SH3 domain located in the cytoplasm. Mutations and deficiencies in peroxisomal protein importing and peroxisome assembly lead to peroxisomal biogenesis disorders, an example of which is Zellweger syndrome. [provided by RefSeq, Oct 2008],</p>