

PX11A Polyclonal Antibody

Cat No: HR1AP11089

For research use only

Overview

Product Name	PX11A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GeneID?Human?	PEX11A PEX11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>peroxisomal biogenesis factor 11 alpha(PEX11A) Homo sapiens This gene is a member of the PEX11 family, which is composed of membrane elongation factors involved in regulation of peroxisome maintenance and proliferation. This gene product interacts with peroxisomal membrane protein 19 and may respond to outside stimuli to increase peroxisome abundance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2012],</p>