

OSBL6 Polyclonal Antibody

Cat No: HR1AP11077

For research use only

Overview

Product Name	OSBL6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	OSBPL6 ORP6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	oxysterol binding protein like 6(OSBPL6) Homo sapiens This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. Transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],