

## S22A4 Polyclonal Antibody

Cat No: HR1AP11072

For research use only

### Overview

Product Name	S22A4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	SLC22A4 ETT OCTN1 UT2H
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 22 member 4(SLC22A4) Homo sapiens Polyspecific organic cation transporters in the liver, kidney, intestine, and other organs are critical for elimination of many endogenous small organic cations as well as a wide array of drugs and environmental toxins. The encoded protein is an organic cation transporter and plasma integral membrane protein containing eleven putative transmembrane domains as well as a nucleotide-binding site motif. Transport by this protein is at least partially ATP-dependent. [provided by RefSeq, Jul 2008].</p>