

NOL7 Polyclonal Antibody

Cat No: HR1AP11068

For research use only

Overview

Product Name	NOL7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	NOL7 C6orf90 NOP27
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nucleolar protein 7(NOL7) Homo sapiens The protein encoded by this gene localizes to the nucleolus, where it maintains nucleolar structure and cell growth rates. The encoded protein also functions as a tumor suppressor and regulator of angiogenesis. The RB tumor suppressor gene recruits transcription factors to this gene and positively regulates its expression. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Dec 2015],