

NOL4 Polyclonal Antibody

Cat No: HR1AP11067

For research use only

Overview

Product Name	NOL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	NOL4 NOLP HRIHFB2255
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tissue specificity:Expressed predominantly in fetal brain, adult brain and testis.,