

NOP14 Polyclonal Antibody

Cat No: HR1AP11066

For research use only

Overview

Product Name	NOP14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	94kDa
GeneID?Human?	NOP14 C4orf9 NOL14 RES4-25
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NOP14 nucleolar protein(NOP14) Homo sapiens This gene encodes a protein that plays a role in pre-18s rRNA processing and small ribosomal subunit assembly. The encoded protein may be involved in the regulation of pancreatic cancer cell proliferation and migration. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014],