

NOL11 Polyclonal Antibody

Cat No: HR1AP11064

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | NOL11 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 79kDa |
| GeneID?Human? | NOL11 L14 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | NOL11 (Nucleolar Protein 11) is a Protein Coding gene. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. Ribosome biogenesis factor. May be required for both optimal rDNA transcription and small subunit (SSU) pre-rRNA processing at sites A, A0, 1 and 2b. |