

NOC4L Polyclonal Antibody

Cat No: HR1AP11061

For research use only

Overview

Product Name	NOC4L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	NOC4L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the CBF/MAK21 family.,