

## SC65 Polyclonal Antibody

Cat No: HR1AP11060

For research use only

### Overview

Product Name	SC65 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	LEPREL4 NOL55 SC65
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>prolyl 3-hydroxylase family member 4 (non-enzymatic)(P3H4) Homo sapiens This nucleolar protein was first characterized because it was an autoantigen in cases on interstitial cystitis. The protein, with a predicted molecular weight of 50 kDa, appears to be localized in the particulate compartment of the interphase nucleolus, with a distribution distinct from that of nucleolar protein B23. During mitosis it is associated with chromosomes. [provided by RefSeq, Jul 2008].</p>