

## STP2 Polyclonal Antibody

Cat No: HR1AP11059

For research use only

### Overview

Product Name	STP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	TNP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:In the elongating spermatids of mammals, the conversion of nucleosomal chromatin to the compact, non-nucleosomal form found in the sperm nucleus is associated with the appearance of a small set of basic chromosomal transition proteins.,similarity:Belongs to the nuclear transition protein 2 family.,tissue specificity:Testis.,