

NRBP2 Polyclonal Antibody

Cat No: HR1AP11058

For research use only

Overview

Product Name	NRBP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	NRBP2 PP9320 TRG16
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	domain:The protein kinase domain is predicted to be catalytically inactive.,sequence caution:Intron retention.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family.,similarity:Contains 1 protein kinase domain.,