

RORB Polyclonal Antibody

Cat No: HR1AP11057

For research use only

Overview

Product Name	RORB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	RORB NR1F2 RZRB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RAR related orphan receptor B(RORB) Homo sapiens The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It is a DNA-binding protein that can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The encoded protein has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, and to help regulate the expression of some genes involved in circadian rhythm. [provided by RefSeq, Feb 2014],</p>