

NU205 Polyclonal Antibody

Cat No: HR1AP11050

For research use only

Overview

Product Name	NU205 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	221kDa
GeneID?Human?	NUP205 C7orf14 KIAA0225
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nucleoporin 205(NUP205) Homo sapiens This gene encodes a nucleoporin, which is a subunit of the nuclear pore complex that functions in active transport of proteins, RNAs and ribonucleoprotein particles between the nucleus and cytoplasm. Mutations in this gene are associated with steroid-resistant nephrotic syndrome. [provided by RefSeq, Jul 2016],