

NF2L1 Polyclonal Antibody

Cat No: HR1AP11044

For research use only

Overview

Product Name	NF2L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	NFE2L1 HBZ17 NRF1 TCF11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nuclear factor, erythroid 2 like 1(NFE2L1) Homo sapiens This gene encodes a protein that is involved in globin gene expression in erythrocytes. Confusion has occurred in bibliographic databases due to the shared symbol of NRF1 for this gene, NFE2L1, and for "nuclear respiratory factor 1" which has an official symbol of NRF1. [provided by RefSeq, Jul 2008],