

NOSIP Polyclonal Antibody

Cat No: HR1AP11040

For research use only

Overview

Product Name	NOSIP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	NOSIP CGI-25
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	nitric oxide synthase interacting protein(NOSIP) Homo sapiens The protein encoded by this gene may modulate the activity and localization of nitric oxide synthase (endothelial and neuronal) and thus nitric oxide production. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Aug 2012],