

## NID2 Polyclonal Antibody

Cat No: HR1AP11038

For research use only

### Overview

Product Name	NID2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	151kDa
GeneID?Human?	NID2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>nidogen 2(NID2) Homo sapiens This gene encodes a member of the nidogen family of basement membrane proteins. This protein is a cell-adhesion protein that binds collagens I and IV and laminin and may be involved in maintaining the structure of the basement membrane.[provided by RefSeq, Jun 2010],</p>