

NETO1 Polyclonal Antibody

Cat No: HR1AP11036

For research use only

Overview

Product Name	NETO1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	NETO1 BTCL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neuropilin and tolloid like 1(NETO1) Homo sapiens This gene encodes a predicted transmembrane protein containing two extracellular CUB domains followed by a low-density lipoprotein class A (LDLa) domain. A similar gene in mice encodes a protein that plays a critical role in spatial learning and memory by regulating the function of synaptic N-methyl-D-aspartic acid receptor complexes in the hippocampus. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2011],</p>