

## GPM6B Polyclonal Antibody

Cat No: HR1AP11031

For research use only

### Overview

Product Name	GPM6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	GPM6B M6B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>glycoprotein M6B(GPM6B) Homo sapiens This gene encodes a membrane glycoprotein that belongs to the proteolipid protein family. Proteolipid protein family members are expressed in most brain regions and are thought to be involved in cellular housekeeping functions such as membrane trafficking and cell-to-cell communication. This protein may also be involved in osteoblast differentiation. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are located on chromosomes Y and 22. [provided by RefSeq, Jan 2016],</p>