

## NAV3 Polyclonal Antibody

Cat No: HR1AP11026

For research use only

### Overview

Product Name	NAV3 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	262kDa
GeneID?Human?	NAV3 KIAA0938 POMFIL1 STEERIN3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neuron navigator 3(NAV3) Homo sapiens This gene belongs to the neuron navigator family and is expressed predominantly in the nervous system. The encoded protein contains coiled-coil domains and a conserved AAA domain characteristic for ATPases associated with a variety of cellular activities. This gene is similar to unc-53, a <i>Caenorhabditis elegans</i> gene involved in axon guidance. Multiple alternatively spliced transcript variants for this gene have been described but only one has had its full-length nature determined. [provided by RefSeq, Jul 2008].</p>