

NMB Polyclonal Antibody

Cat No: HR1AP11025

For research use only

Overview

Product Name	NMB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	NMB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neuromedin B(NMB) Homo sapiens This gene encodes a member of the bombesin-like family of neuropeptides, which negatively regulate eating behavior. The encoded protein may regulate colonic smooth muscle contraction through binding to its cognate receptor, the neuromedin B receptor (NMBR). Polymorphisms of this gene may be associated with hunger, weight gain and obesity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2015],</p>