

NLGN3 Polyclonal Antibody

Cat No: HR1AP11024

For research use only

Overview

Product Name	NLGN3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	93kDa
GeneID?Human?	NLGN3 KIAA1480 NL3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neuregulin 3(NLGN3) Homo sapiens This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. Mutations in this gene may be associated with autism and Asperger syndrome. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Oct 2009],</p>