

NLGN1 Polyclonal Antibody

Cat No: HR1AP11022

For research use only

Overview

Product Name	NLGN1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	NLGN1 KIAA1070
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neuroligin 1(NLGN1) Homo sapiens This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. [provided by RefSeq, Jul 2008],</p>