

NGDN Polyclonal Antibody

Cat No: HR1AP11021

For research use only

Overview

Product Name	NGDN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	NGDN C14orf120
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	neuroguidin(NGDN) Homo sapiens Neuroguidin is an EIF4E (MIM 133440)-binding protein that interacts with CPEB (MIM 607342) and functions as a translational regulatory protein during development of the vertebrate nervous system (Jung et al., 2006 [PubMed 16705177]).[supplied by OMIM, Mar 2008],