

NDF4 Polyclonal Antibody

Cat No: HR1AP11017

For research use only

Overview

Product Name	NDF4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	NEUROD4 ATH3 ATOH3 BHLHA4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Appears to mediate neuronal differentiation.,PTM:Serine or threonine phosphorylation within the basic region may regulate neurogenic activity.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,