

NDF6 Polyclonal Antibody

Cat No: HR1AP11018

For research use only

Overview

Product Name	NDF6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	NEUROD6 ATOH2 BHLHA2 My051
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	neuronal differentiation 6(NEUROD6) Homo sapiens This gene is a member of the NEUROD family of basic helix-loop-helix transcription factors. The encoded protein may be involved in the development and differentiation of the nervous system. [provided by RefSeq, Nov 2012],