

NRSN2 Polyclonal Antibody

Cat No: HR1AP11007

For research use only

Overview

Product Name	NRSN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	NRSN2 C20orf98
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:It is uncertain whether Met-1 or Met-2 is the initiator.,function:May play a role in maintenance and/or transport of vesicles.,similarity:Belongs to the VMP family.,