

CHL1 Polyclonal Antibody

Cat No: HR1AP11006

For research use only

Overview

Product Name	CHL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	132kDa
GeneID?Human?	CHL1 CALL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cell adhesion molecule L1 like(CHL1) Homo sapiens The protein encoded by this gene is a member of the L1 gene family of neural cell adhesion molecules. It is a neural recognition molecule that may be involved in signal transduction pathways. The deletion of one copy of this gene may be responsible for mental defects in patients with 3p- syndrome. This protein may also play a role in the growth of certain cancers. Alternate splicing results in both coding and non-coding variants. [provided by RefSeq, Nov 2011],</p>