

# UNC5C Polyclonal Antibody

Cat No: HR1AP11003

For research use only

## Overview

Product Name	UNC5C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	UNC5C UNC5H3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>unc-5 netrin receptor C(UNC5C) Homo sapiens This gene product belongs to the UNC-5 family of netrin receptors. Netrins are secreted proteins that direct axon extension and cell migration during neural development. They are bifunctional proteins that act as attractants for some cell types and as repellents for others, and these opposite actions are thought to be mediated by two classes of receptors. The UNC-5 family of receptors mediate the repellent response to netrin; they are transmembrane proteins containing 2 immunoglobulin (Ig)-like domains and 2 type I thrombospondin motifs in the extracellular region. [provided by RefSeq, Jul 2008].</p>