

UNC5B Polyclonal Antibody

Cat No: HR1AP11002

For research use only

Overview

Product Name	UNC5B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	103kDa
GeneID?Human?	UNC5B P53RDL1 UNC5H2 UNQ1883/PRO4326
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>unc-5 netrin receptor B(UNC5B) Homo sapiens This gene encodes a member of the netrin family of receptors. This particular protein mediates the repulsive effect of netrin-1 and is a vascular netrin receptor. This encoded protein is also in a group of proteins called dependence receptors (DpRs) which are involved in pro- and anti-apoptotic processes. Many DpRs are involved in embryogenesis and in cancer progression. Two alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Oct 2011],</p>